

i-STAT SYSTEM: RECEIPT OF NEW CARTRIDGES & MONTHLY CHECK

- Test one cartridge from each lot in every shipment upon receipt.
- 4 Window Temperature Strip → If the windows "3" or "4" are red, do not use the cartridges or controls.
- WarmMark® Temp Tag → If window is red, indicates temp above range, do not use cartridges or controls.
- Test a single cartridge from the refrigerator monthly.
- Control ampule must be at room temperature for at least 30 minutes prior to using.
- Ensure Value Assignment Sheet matches the control lot and expiration date to determine acceptable range.

MONTH: _____ YEAR: _____ PRACTICE LOCATION: _____

RECEIPT OF NEW CREATININE CARTRIDGES (Unit: mg/dL)

RECEIVED DATE:				
TESTED BY:				
CLEW:				
CARTRIDGE TYPE:	Crea	Crea	Crea	Crea
CARTRIDGE LOT #:				
CARTRIDGE EXPIRATION DATE:				
QUANTITY:				
CONTROL TYPE & LEVEL	i-STAT Level 1	i-STAT Level 1	i-STAT Level 1	i-STAT Level 1
CONTROL LOT #:				
CONTROL EXPIRATION DATE:				
TEMPERATURE STRIP COLOR:				
ACCEPTABLE RANGE: (from VAS)				
RESULTS:				
RESULTS IN RANGE?				
RESULTS FOR REPEAT TEST: (if applicable)				
REPEAT RESULTS IN RANGE?				

CREATININE CARTRIDGE MONTHLY CHECK (Unit: mg/dL)

DATE	TESTED BY	CLEW	CARTRIDGE LOT #	CONTROL LEVEL	CONTROL LOT #	RANGE	RESULTS	IN RANGE?
				i-STAT Level 1				